

## manual vacuum aspirator

### Technical Specifications

Aspirate Container Capacity: 240 ml (excluding cap)

Occluded suction: -40 kPa within 10s @ less than 2Hz  
pump rate

Internal Diameter of Catheter mount: 17mm

Overflow safeguard: Aspirator functions normally with  
automatic aspirate ejection.

Operating Principle: Rolling seal diaphragm

Design Life: 5 years with annual checkups

Operating Temperature: -20 to 50°C

Dimensions: Net 165 x 165 x 87 mm

Gross 212 x 102 x 233 mm

Weight: Net 400g

Gross 650g

### **Contact us**

Steel yao

Meds Medical Equipment Co.,LTD

Suite 808,Taikun International Plaza,

Jiaxing,Zhejiang,China.

Tel:+86-57383768915

MP:+86-15857396317

Fax:+86-57383667955

Email:[steel@medinstru.com](mailto:steel@medinstru.com)

website: [www.medinstru.com](http://www.medinstru.com)